

Figure 1 CMS Logo



Measure Authoring Tool Release Notes

Version 4.4.2

May 10, 2016

Table of Contents

ENHANCEMENTS	3
1.1 SECURITY UPDATES	3
1.2 ALLERGY UNIT (AU) AND BIOEQUIVALENT ALLERGEN UNIT (BAU) ADDED	3
RESOLVED DEFECTS	4
2.1 PACKAGING ERROR WITH DEEPLY NESTED CLAUSES RESOLVED	4
2.2 TAG UPDATED WITH MORE DESCRIPTIVE TEXT WHEN RESETTING PASSWORDS USING ASSISTIVE DEVICES	4

Enhancements

Measure Authoring Tool version 4.4.2 enhancements are described here.


1.1 Security Updates

Annually, the MAT undergoes a security audit to ensure the MAT security features continue to evolve as security recommendations and practices evolve. As a result there are a few minor updates implemented with this version of the MAT to incorporate additional security precautions. There are no changes to the MAT features and functionality with this security update.

1.2 Allergy Unit (AU) and Bioequivalent Allergen Unit (BAU) Added

The dosing units Allergy Unit (AU) and Bioequivalent Allergen Unit (BAU) are added to the available units in the Clause Workspace when adding attributes.

Figure 2 AU and BAU Added



Resolved Defects

2.1 Packaging Error with Deeply Nested Clauses Resolved

A packaging error was reported with the use of deeply nested clauses (more than 6 levels). Users were prevented from completing a measure package, despite the included clauses and population logic successfully passing validation(s). This error is resolved.

2.2 Tag Updated with More Descriptive Text When Resetting Passwords Using Assistive Devices

To ensure individuals requiring the use of assistive devices are able to successfully reset a password using the password reset option in the MAT Account tab, the tag is updated. The update includes more descriptive verbiage and reordering of alternative text so a user using assistive technology may successfully enter and confirm a new password.