## Measure Authoring Tool 6.10 QDM and FHIR

MAT Maintenance Release September 29, 2021





## Included in the MAT 6.10 Maintenance Release:

- The Timeout in the MAT has increased from 15 to 30 minutes
  - •The Inactivity Warning has increased from 10 to 25 minutes
- For FHIR Measures and Libraries, the CQL to ELM Translator has been upgraded to version 1.5.4
  - •QDM Measures and Libraries still use version 1.5.3
- To address potential security vulnerabilities:
  - •MAT was upgraded to Tomcat 10
  - •MAT microservices were upgraded to Java 15

## **Contact Us**

The MAT can be accessed from the MAT public web site <a href="here">here</a>.

Detailed release notes are located <u>here</u> on the MAT public web site, under the Training and Resources page.

The MAT support team is available to answer questions and to provide user support.

Please continue to contact MAT Support through the Bonnie and MAT Issue Tracker here: <a href="https://oncprojectracking.healthit.gov/support/projects/BONNIEMAT/summary">https://oncprojectracking.healthit.gov/support/projects/BONNIEMAT/summary</a>.