



Measure Authoring Tool 6.10 QDM and FHIR

*MAT Maintenance Release
November 10, 2021*

Included in the MAT 6.10 Maintenance Release:

- The Measure Export FHIR Bundles will now include entries for nested CQL Libraries. Previously, only CQL libraries at the top level were included in the FHIR Bundle.
- The Human Readable export will now only include the Measure Observation section in the header when appropriate. Previously, the Measure Observation section was displayed regardless of measure type.

Contact Us

The MAT can be accessed from the MAT public web site [here](#).

Detailed release notes are located [here](#) on the MAT public web site, under the Training and Resources page.

The MAT support team is available to answer questions and to provide user support.

Please continue to contact MAT Support through the Bonnie and MAT Issue Tracker here: <https://oncprojecttracking.healthit.gov/support/projects/BONNIEMAT/summary>.