

MADiE 1.1.12 Release Notes

The following updates are included in the MADiE 1.1.12 release:

Category	Feature Description
Test Case Creation	QI-Core Measures: MADiE has been updated to provide the ability to delete all the test cases associated with the measure on the test case list screen by pressing the Delete All button. The user will be asked to confirm this action before removing all the test cases from the measure.
Test Case Export	QI-Core Measures: MADiE has been updated to enforce a unique test case group and title combination across the measure when importing test cases. Test cases without a unique group and title combination will not be imported.
Test Case Export	QI-Core Measures: MADiE has been updated to remove leading and trailing spaces and spaces within test case title and group values when the test case filename is created during test case export. BONNIEMAT-1550
UI Enhancements	<p>QI-Core Measures: Several improvements have been made to the MADiE User Interface on the edit test case screen</p> <ul style="list-style-type: none"> The Highlighting tab and the CQL code now have the same height and allow independent scrolling which lets the user compare any part of the CQL to any part of the Highlighting content in a side-by-side manner. BONNIEMAT-1495 The Validation Errors panel has been expanded to have the same height as the other panels for a more consistent user experience. Breadcrumbs have been added to the top of the edit test case screen to improve navigation. The breadcrumbs will have a link back to the test case list page and it will also display the group and title of the test case in which you are currently working.
Measure Creation	The "Program Use Context" field has been removed from the Details tab for both measures and libraries. BONNIEMAT-1533
Measure Creation	MADiE has been updated to reduce the character limit to 64 characters instead of 255 characters. The Library Name field on the Create/Edit Measure screens and Create/Edit Library screens have been updated.
Measure Export	QI-Core Measures: MADiE has been updated to use the latest LiquidBase Templates based on QMIG profiles. More information can be found here . BONNIEMAT-1556
Resolved Known Issues	<p>Date Added: 4/12/2023 BONNIEMAT-1429 BONNIEMAT-1452</p> <p>QI-Core Measures: Measure export may fail when text fields include some symbols. Less than (<) and ampersand (&) symbols have been identified as causing export to fail if they are included, but additional symbols may have this effect as well.</p> <p>BONNIEMAT-1551</p> <p>Previously, paging through many pages of measures on the All Measures tab would eventually result in an error. This has been fixed and now you can page through the full list of measures.</p>

MADiE 1.1.12 Release Notes

Date Added: 8/23/2023 [BONNIEMAT-1584](#)

The MADiE 1.1.10 release included updates to HAPI which resulted in some test cases becoming invalid due to having validation errors. The test case status that displays on the test case list page is not currently updated until after each test case is updated and saved. MADiE users may find their test case list displaying all test cases as valid, when one or more test cases may no longer be valid and invalid test cases are not included in test case exports.

MADiE will run validations on all draft measure test cases when a release includes data or HAPI updates. This will ensure the test case list page displays the updated status for each test case after the release is complete.