

MADiE 1.3.0 Release Notes

The following updates are included in the MADiE 1.3.0 release:

Category	Feature Description
CQL Libraries	MADiE now allows users to share edit access to stand-alone CQL Libraries. To share edit access of a library with another user, fill out the MADiE Measure or Library Sharing Request Form and email it to the help desk at semanticbits-madie-help@icf.com .
Measure Export	QDM Measures: MADiE has been updated to allow users to export QDM measures from MADiE. Please note that measures transferred from MAT to MADiE need to be reviewed in MADiE and the CQL updated (can just add a space to the CQL) after the transfer and prior to working with the measure (including working with test cases or exporting the measure). Until this step is complete you may see a CQL error notification on the Test Cases screen and/or during Measure Export.
Test Case Creation	QDM Measures: Users can now choose to include Supplemental Data Elements (SDE) in their test case calculations and highlighting.
Measure Export	MADiE has been updated to use LiquidBase Templates 0.4.1 based on QMIG profiles. More information can be found here .
Resolved Known Issues	Date Added: 3/15/24 QDM Measures: When a measure is retransferred from the MAT to MADiE, if the measure had test cases added or imported in MADiE prior to retransferring the measure, the test case expected values in MADiE are deleted during the retransfer. BONNIEMAT-1772 QDM Measures: MADiE is filtering out DRCs that are not part of Source Data Criteria when building out the CQM Measure's list of value sets. This has been fixed and the DRCs that are not part of Source Data Criteria are available now. BONNIEMAT-1770