

MADiE 2.0.0 Release Notes

The following updates are included in the MADiE 2.0.0 release:

| Category | Feature Description |
|-----------------------|--|
| UI Enhancements | When a CMS Web Application Firewall (WAF) error blocks an action from being performed, a popup message will now display. This message will notify the user they have encountered a CMS WAF issue and instruct the user to create a MADiE IT ticket and, if available, will display the support ID to include in the ticket. |
| CQL-ELM Translator | <p>MADiE has been updated to support updates being made to the CQL to ELM translator version for QDM or for QI-Core separately.</p> <p>MADiE now displays the version of the CQL to ELM translator that is used in the measure.</p> <ul style="list-style-type: none"> Measures in Draft status – the "Currently using CQL to ELM Translator Version" field was added to display the current translator version Measures in Versioned status – the "Versioned with CQL to ELM Translator Version" field was added to display the version of the translator that was used when versioning the measure |
| Test Case Execution | QDM Measures: The Manifest selected for value set expansion will be used for testing. Value set versions specified in the measure CQL cannot be used and the value set version will be removed from the CQL the next time the user saves the measure from any tab. After saving, a message will display notifying the user of this update to the measure CQL. |
| Measure Export | QDM Measures: The list of values in the Terminology section of the human readable export for QDM measures is now correctly sorted (BONNIEMAT-1868) BONNIEMAT-1946) |
| Measure Creation | QI-Core Measures: The stratification drop down has been updated to include <i>only</i> the populations selected in the population criteria. |
| Resolved Known Issues | <p>Date Added: 6/20/2024</p> <p>QDM Measures: Users may experience unexpected behavior with the test case “Delete All” button.</p> <ul style="list-style-type: none"> The “Delete All” button will remain enabled after users delete all test cases. The “Delete All” button will remain disabled after importing test cases. <p>Date Added: 6/20/2024</p> <p>QDM Measures: After importing test cases, an error will display on the test case list page: “An error occurred, please try again. If the error persists, please contact the help desk.”</p> |

MADiE 2.0.0 Release Notes

| Category | Feature Description |
|--------------------------|--|
| Resolved Known Issues | <p>Date Added: 6/20/2024 QDM Measures: Initially after users change the selected Manifest for test case execution, the test case results may display incorrect information.</p> <p>Date Added: 6/20/2024 QDM Measures: When users edit a test case, select “Run Test” and then select “Discard Changes” the Expected/Actual tab may display incorrect results as the results will revert to the last saved expected and actual values. The results displayed on the test case list page are not impacted.</p> <p>Date Added: 6/20/2024 (BONNIEMAT-1894, BONNIEMAT-1933) QDM Measures: When updating Stratification definitions (the narrative) after transferring a measure from MAT to MADiE, users may receive the following error message. "QDM group stratifications cannot be associated." Work around: Delete the Stratification groups and rebuild in MADiE.</p> |