

MADiE 3.0.5 Release Notes

The following updates are included in the MADiE 3.0.5 release:

| Category | Feature Description |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Measure Creation | <p>CMS IDs are now displayed and exported as four-digit values, with leading zeros added when needed (for example, “2” appears as “0002”).</p> <p>Note: While all measures show the padded CMS ID in the UI, only export packages generated from draft measures, or versions created after this release, include the four-digit format. Exports created from measures versioned before this release remain unchanged. (MADIE-2364)</p> |
| Test Case Creation | <p>QI-Core Measures: A QI-Core Patient will now be automatically added by default to the Added tab and the test case JSON when creating a new test case.</p> |
| CQL-ELM Translator | <p>MADiE has been updated to use CQL-ELM Translator 4.8.0 for both QDM and QI-Core. More information can be found here.</p> |
| QI-Core Liquid Template | <p>MADiE has been updated to use Liquid Templates 0.5.5. More information can be found here.</p> |
| HAPI FHIR | <p>MADiE has been updated to use HAPI FHIR version 8.10.0. More information can be found here.</p> |
| Administrator Functionality | <p>MADiE Administrator functionality has been updated.</p> |